

Gregory Seat Technology



Standard Deluxe

Gregory's Deluxe Seat uses moulded foam with anti-bacteriaside in it to prevent the spread of infection in the seat foam. The contour of the seat is dished to form to the shape of the user with a waterfall at the front which alleviates pressure at the front of the thigh.

Why it works:

- Back of the seat is contoured so you sit into the seat.
- Pommel at the front to direct the legs in the correct position.
- Waterfall at the front of the seat to reduce pressure under the thighs.
- Anti-bacterial moulded foam

